



FACILITIES

Emtech Laboratories
7745 Garland Circle
Roanoke, VA 24019

MAILING

Emtech Laboratories
P.O. Box 12900
Roanoke, VA 24022

PHONE

Local: (540) 265-9156
Toll Free: (800) 336-5719
Fax: (540) 265-9164

M.C.T.

M.C.T. earpieces are great for high frequency losses and other losses not greater than 60db at 1kHz.

All but eliminates the "voice in the barrel" (occlusion effect) problem.
Requires less venting, thereby reducing the risk of acoustic feedback.

Reduces contact with the walls of the ear canal for increased comfort.

If canal retention and/or patient dexterity are not a problem most M.C.T.'s will be made as a canal only, producing a more cosmetically appealing earpiece.

The canal length in M.C.T. impressions must be 15mm to 20mm.

