



FACILITIES
Emtech Laboratories
7745 Garland Circle
Roanoke, VA 24019

MAILING
Emtech Laboratories
P.O. Box 12900
Roanoke, VA 24022

PHONE
Local: (540) 265-9156
Toll Free: (800) 336-5719
Fax: (540) 265-9164

Semi Skeleton

This style is basically a skeleton with the upper portion of the concha ring removed. In most cases it offers as much retention as the skeleton and is lighter. This is an excellent selection when the ear is not undercut in the upper concha and full retention can be obtained in the helix and anti-tragus.

