



FACILITIES  
Emtech Laboratories  
7745 Garland Circle  
Roanoke, VA 24019

MAILING  
Emtech Laboratories  
P.O. Box 12900  
Roanoke, VA 24022

PHONE  
Local: (540) 265-9156  
Toll Free: (800) 336-5719  
Fax: (540) 265-9164

## Skeleton

This style offers excellent retention, yet is very lightweight. It is completely cut out, leaving only the outer ring it offers good concealment as well as comfort and comes in the material of your choice. Skeleton style is also used for the communication earpieces.

