



FACILITIES
Emtech Laboratories
7745 Garland Circle
Roanoke, VA 24019

MAILING
Emtech Laboratories
P.O. Box 12900
Roanoke, VA 24022

PHONE
Local: (540) 265-9156
Toll Free: (800) 336-5719
Fax: (540) 265-9164

TV Communication Ear Pieces/ Molds

There are 3 different styles of TV Communication Molds available.

- 1- Custom fitted with standard tubing is an ear piece made of Lucite in a skeleton style, non-occluding with a straight piece of tubing and female adapter attached.
- 2- Custom fitted with coiled tubing is the same as number one except it has a coiled tubing, allowing more freedom of movement.
- 3- Universal one-size-fits-all ear tips with straight tubing, female adapter and Comply foam plug.

