



FACILITIES
Emtech Laboratories
7745 Garland Circle
Roanoke, VA 24019

MAILING
Emtech Laboratories
P.O. Box 12900
Roanoke, VA 24022

PHONE
Local: (540) 265-9156
Toll Free: (800) 336-5719
Fax: (540) 265-9164

TechSeal

The TechSeal permits boosting high frequency sound without causing aggravating feedback experienced with other earpieces. Its versatility is unmatched, covering both severe and profound losses, as well as mild and moderate losses and acoustic modifications.

The canal of the earpiece is hollow. Combine that with quality in a mold that is so soft and pliable that you may not even realize it's in the ear. An Innovative Tube lock holds the tubing securely without cement. (Very easy to retube with no mess and no waiting for glue to dry.) Standard canal with half lock makes it very light and comfortable, and will not work its way out of the ear.

Other styles available upon request.

